

## Clarified STAR Protocol: Antimicrobial susceptibility testing to evaluate minimum inhibitory concentration values of clinically relevant antibiotics




Figure 1. MIC Assay Flow Through

### Before you begin:

- I. Determine the optimal MIC testing range for your pathogen using the EUCAST MIC distribution repository (<https://mic.eucast.org/search/>).
  - i. Select a continuous range of ten 2-fold serial dilutions that would encompass the clinical breakpoints used to categorize the pathogens as susceptible, intermediate, or resistant.
  - ii. Calculate 2x the highest antibiotic concentration based on this range for the maximum concentration for your serial dilutions.
- II. Prepare your stock solutions and culture media in the correct solvents (see the list on the fridge for reference). Store at 4°C for short term storage or -20°C for long term storage.
  - i. Sterile LB (with/without antibiotics) ~ 15-50 mL
  - ii. LB + 2x highest antibiotic concentration ~ 2-5 mL

### Protocol:

Timing: 2 Days

1. Streak out culture of choice for isolation on an LB + Antibiotic plate and grow at 37°C for 18 hrs.
2. Inoculate 3 mL sterile LB without antibiotics with a single colony from your isolation plate.
  - a. If you have a strain that has been selected for using antibiotics, continue to include it in your culture medium.
  - b. Grow overnight for 18 hours at 37°C and 220 rpm.

Timing: 2-3 hours

3. Using your prepared LB + 2x antibiotic stock, add 100 µL of your stock to wells A<sub>1</sub>-H<sub>1</sub> of a clear 96 well plate.
4. Add 50 µL of LB, or appropriate media, to each well in columns 2-12.
  - a. Add an additional 50 µL of LB to column 12 for your negative control (media only).

### Serial Dilutions:

5. Transfer 50  $\mu\text{L}$  from column 1 into column 2, pipetting up and down 3 times to mix.
  - a. Repeat the transfer process in columns 3-10.
  - b. Discard the additional 50  $\mu\text{L}$  of culture from column 10.




Figure 2. [Serial Dilutions of Antibiotics for MIC Assay](#)



### Culture Dilutions and Final Plating:

6. Add 100  $\mu\text{L}$  of your 18-hour overnight cultures to a different clean clear 96 well plate. Include a LB blank. Measure the OD600.
  - a. The spectrophotometer measures both alive and dead cells in a broth culture. Because of this a standard line needed to be created to calculate the CFU/mL at a given absorbance value.
  - b. Create 6-8 dilutions of your culture and measure their absorbances. Plot these values against CFU (cells/mL)  $\times 10^6$ . With the slope of the line, you can determine the CFU/mL at the experimental OD.
  - c. For a full explanation on determining CFU/mL, see this link here → [1.15: Determination of Bacterial Numbers - Biology LibreTexts](#)
7. Dilute your overnight cultures to  $10^6$  CFU/mL in LB (2x the desired bacterial concentration).
8. Add 50  $\mu\text{L}$  of your diluted cultures to all wells on the plate except for those in column 12.
9. Incubate plates at 37°C without shaking for 20 hours.

## Determine Experimental MIC:


Timing: 20 min

10. Score the plate based on
  - a. Presence/absence of turbidity in test wells
  - b. Confirm growth in bacteria only positive control wells
  - c. Confirm no growth in media only negative control wells
11. Interpret the MIC values with respect to the clinical breakpoints determined through EUCAST.



Kanamycin / *Escherichia coli*  
 International MIC distribution - Reference database 2024-08-31  
**Based on aggregated distributions**

MIC distributions include collated data from multiple sources, geographical areas and time periods and can never be used to infer rates of resistance



MIC  
 Epidemiological cut-off (ECOFF): (16) mg/L  
 Wildtype (WT) organisms: ≤ 16 mg/L

Confidence interval: 4 - 32  
 3860 observations (4 data sources)